

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3

Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz DL156BH - Ultrasound bath DL 156 BH Order code: **2007.7181**



Information about product price on demand

Parameters

Quantitative unit ks

Detailed description

Description

- Cleaning of sensitive materials
- Mixing of samples
- Homogenizing of samples for analysis
- Fast degassing of solvents for HPLC

Features

- Power control from 20 to 100 %
- Heating from 20 to 80 °C, thermostatically adjustable, LED display
- Optical LED signal in case of excess temperature
- Automatic short-time extra ultrasound during heating up period to avoid delay in boiling switch-off possible
- Timer for countdown or continuous operation
- DEGAS quick degassing
- SweepTec® permanent sound field oscillation ensures even and gentle cleaning
- Automatical safety shut-down 12 hours after last keypress
- Marking of filling level for safe dosage
- Drain valve for simple emptying, one-piece welded drain bow
- · Protection class IP 33 at highest possible level through liquid-tight mains cable and housing front without openings

For more information about this product please download product sheet under "Files for download".

Do you want to print the page now? Hide some parts of page before printing



Print partially

Print the page as you see it now



Print complete page

(this box will not be printed)

Edit before printing

Hover your mouse over some part of page, a red rectangle will appear. Click it to hide the selected content of page.



Print!

