



HELAGO-CZ, s.r.o.  
Commercial Register maintained by the Regional Court in Hradec  
Králové  
Section C, File 17879  
Kladská 1082  
500 03 Hradec Králové 3  
Company ID: 25 96 39 61, VAT: CZ 25963961  
Phone: 495 220 229, 495 220 394  
Fax: 495 220 154  
GSM gate: 602 123 096  
E-mail: info@helago-cz.cz  
Web: http://www.helago-cz.cz

**AR1000 - Superior Female  
Bedford Doll**  
Order code: **4101.AR1000**



Information about product price on demand

Parameters

Nursing simulators - filter

Dospělý

Quantitative unit

ks

Part of the Bedford Nurse Training Doll range comprising of seven models available in both Caucasian and Black Skin tones.

- Bedford Nurse Training Dolls are the only models for with a soft and palpable skin which can be washed
- Bedford Nurse Training Dolls are the only models with continuous skin in elbow and knee joint

This model is designed for Nurse Training and offers the following features:

- Life Size measuring 164cms tall
- Full range of movement through neck, shoulder, elbow wrist, hip, knee and ankle joints

- Eyes - for washing and administering eye drops
- Ears - for aural hygiene
- Vagina - for vaginal douching
- Rectum - for enemas
- Soft feet with separate flexible toes for foot hygiene
- Washable, nylon, rooted hair for hair care
- Mouth and nose - for nasal and gastric feeds, aspiration and stomach irrigation
- Stomach reservoir - connected to mouth, nose and gastrostomy
- Urethral reservoir - connected to urethra and supra-pubic orifice
- Rectal Reservoir - connected to anus and colostomy
- Interchangeable male and female genitalia
- Improved male genitalia now has the urethral canal formed in such a realistic way that the penis must be raised to the correct position to allow insertion of catheter.
- Tracheostomy
- Replaceable hypodermic injection pad in upper left arm
- Replaceable intramuscular injection pads in upper left leg and upper outer quadrant of left buttock
- Articulated waist allows the body to twist, bend forwards and backwards, increasing flexibility
- The body separates at the waist which allows access to the internal reservoirs and the arm and leg connectors
- Moulded colostomy stoma for greater realism

This model allows the practice of the following clinical procedures:

- Bed positioning
- Bandaging - both wet and dry
- Washing
- Catheterisation/Bladder irrigation
- Intramuscular injections
- Enemas
- Naso-gastric feeding
- Tracheostomy Care

#### **Optional Accessories**

- AR1000-B - Superior Female Bedford Doll, Black
- AR12 - Above Elbow Amputation Stump
- AR13 - Below Elbow Amputation Stump
- AR14 - Above Knee Amputation Stump
- AR15 - Below Knee Amputation Stump