



HELAGO-CZ, s.r.o.  
Commercial Register maintained by the Regional Court in Hradec Králové  
Section C, File 17879  
Kladská 1082  
500 03 Hradec Králové 3  
Company ID: 25 96 39 61, VAT: CZ 25963961  
Phone: 495 220 229, 495 220 394  
Fax: 495 220 154  
GSM gate: 602 123 096  
E-mail: info@helago-cz.cz  
Web: http://www.helago-cz.cz

## AC8 - Adult Intubation Trainer

Order code: **4101.AC8**



Information about product price on demand

Parameters

Quantitative unit

ks

- This trainer has been developed to learn how to perform intubation under lifelike conditions
- If during intubation, too much pressure is applied through a laryngoscope to the upper incisors, an audio signal is heard
- The correct location of the tube, inserted after intubation, may be confirmed by audible respiratory sounds through a stethoscope (not supplied) lateral to the left or right wall of the thorax
- Includes separate model of the larynx for identification of anatomical structures
- Supplied in a carrying case