


HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394 Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

## 50124 - Vermiform Appendix -Retrocaecal Order code: 4104.50124



Cena bez DPH Price with VAT

Parameters Surgery - filter Quantitative unit 162,80 Eui 196,99 Eur

> Basic surgery ks

For practising the laparoscopic removal of the vermiform appendix. There are 3 variations of this model which present the appendix in the following positions:

- normal (as found in 32% of patients)
- post-ileal (as found in 1% of patients)
- retrocaecal (as found in 64% of patients)

Each model has the additional feature of a simulated peptic ulcer, for closure using peritoneum and mesentry.

## Skills

- Excision and division of peritoneum
- Identification and exposure of appendix
- Mobilisation and division of vessels
- Removal of appendix
- Inspection of stump of appendix
- Repair and closure of peptic ulcer

## Features

- Highly realistic model of appendix, caecum, and ileum stump presented anatomically
- The variations available in the range of appendices provide increased levels of difficulty
- Convincing fluid filled vessels, caecum contents, and tissues encourage correct handling of tissue
- Vessel incorporates a luer lock to which the Mock Blood Giving Set may be attached to provide fluid flow for the appendicular artery. Use Mock Blood -Arterial (250 ml) - **60653**
- Ideal for training course requirements; the model is designed for swift attachment to and removal from the Soft Tissue Retaining Set