



Cena bez DPH

1.529,00 Eur

Price with VAT

1.850,09 Eur

Parameters

Quantitative unit

ks

- Detailed description
- Rotation and flame sterilization happen automatically. No need for cooling down of inoculating loops in agar.
- When needed - at the push of a button - all 5 inoculating loops can be flame sterilized in one full rotational cycle. The time for flame sterilization can continuously be set from 2 to 10 seconds. Suited for right- and left-handed use.

- Can be used with all standard gas safety-burners (adapters are optionally available).
- Flame sterilization time can be minimized to save gas and to increase lifetime of inoculating loops
- No accidental usage of hot loops. Maintenance-free switches of the flame sterilization carousel are covered inside the device, no malfunction due to dirt or stray light.
- We recommend schuett manufactured inoculating loops made of special stainless steel wire or platinum-iridium wire (90/10). Length 60 mm, diameter of wire 0,5 mm.

Do you want to print the page now?

Hide some parts of page before printing



Print partially

Print the page as you see it now



Print complete page

(this box will not be printed)

Edit before printing

Hover your mouse over some part of page, a red rectangle will appear. Click it to hide the selected content of page.



Print!



Show hidden