



HELAGO-CZ, s.r.o.  
Commercial Register maintained by the Regional Court in Hradec  
Králové  
Section C, File 17879  
Kladská 1082  
500 03 Hradec Králové 3  
Company ID: 25 96 39 61, VAT: CZ 25963961  
Phone: 495 220 229, 495 220 394  
Fax: 495 220 154  
GSM gate: 602 123 096  
E-mail: info@helago-cz.cz  
Web: http://www.helago-cz.cz

**LF04023 - KERi Auscultation  
Manikin**  
Order code: **4107.LF04023**



Information about product price on demand

Parameters

Nursing simulators - filter

Quantitative unit

Dospělý

ks

For auscultation training, the instructor selects from a menu of heart and lung conditions by wireless remote control. The manikin presents itself as a real patient without visible auscultation sites. The student must palpate to identify correct auscultation locations, and will hear different heart and lung sounds as the SmartScope™ is moved to different locations on the manikin. A diagnosis of the condition selected by the instructor can be made by comparing the variations in sounds occurring at different sites - just like a real patient! Lung sounds can be detected at five anterior and 12 posterior locations and students can practice auscultation at six anterior heart sites. The remote control does not have to be pointed directly at the manikin or SmartScope™ to operate. One remote control will operate multiple sets of SmartScopes™ and manikins simultaneously. Great for group instruction. Listen to the sounds by using either the single- or dual-user

headpieces on the SmartScope™ or connecting to an amplified speaker (not included, see SB20146U). Range on the unit is up to 100 ft. of remote access. Each trainer includes manikin, one remote control with LCD display, and one SmartScope™ with single- and dual-user headpiece. Requires two “AA” and two “AA” batteries (included).

**Has all the same features as the Complete KERi™ Manikin (LF04020U) and:**

Heart Conditions:

- Normal
- Aortic Regurgitation
- Mitral Stenosis
- Pulmonary Stenosis
- Holosystolic
- Mid Systolic
- S3 Gallup
- S4 Gallup
- Systolic Click
- Atrial Septal Defect
- PDA
- VSD

Lung Conditions:

- Normal
- Wheeze
- Mono Wheeze
- Fine Crackle
- Rhonchi
- Stridor
- Bronchial
- Coarse Crackle
- Cavernous
- Pulmonary Edema
- Friction Rub
- Pectoriloquy