



HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional Court in Hradec Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961
Phone: 495 220 229, 495 220 394
Fax: 495 220 154
GSM gate: 602 123 096
E-mail: info@helago-cz.cz
Web: http://www.helago-cz.cz

Decubitus Treatment Simulator

Order code: **4109.M156**



Cena bez DPH

2.178,00 Eur

Price with VAT

2.635,38 Eur

Parameters

Quantitative unit

ks

(4 types in 1 set)

Directors: Takaoka station south clinic president Kunio Tsukada, Kanazawa University premedical course Hiromi Sanada

Decubitus series describing particularity of each stage from I-IV with reality.

By fastening to the human body or a model doll, practical training can be done on treatment of Decubitus. There are 4 types (stages I-IV) in 1 set. Colors have been added to express the different features of each stage. They have been done from the inside, so the colors will not drop even if chemicals are used or when washed. The colors of the part infected by Decubitus will not change even after the drill. The structure of damage in the epiderm and the soft touch of each stage have been reproduced realistically, so one can practice treatment in the various stages. It comes with a fastening belt making it possible to attach to a human body.

Stage I. The epiderm has not yet been torn but a rash can be observed.

Stage II. Abrasion, bleeding, and effusion can be observed.

Stage III. Healing process of stage IV where the fascia and deeper had been affected.

Stage IV. There exists a hole widely separating the fascia and periosteum from the hypodermal organ.