

**Wax Form (class III)**  
Order code: **4112.345149**



Cena bez DPH

50,40 Eur

Price with VAT

60,98 Eur

Parameters

Quantitative unit

ks

- Wax form based on the classification of Angle malocclusion

**Examples of Practicing with the Typodont Practice Model:**

- Attach the ORT various Wax Form to the Typodont Ocluder by applying molten wax with an instrument. (use the included wax block)
- Remove the jaws with the black base
- Dip the jaws into 55°C~60°C warm water
- Lift the jaws from the warm water and dip into cool water. Observe the result