

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394 Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

## Eppendorf Reference ® 2, variable 10 - 100 μl Order code: 1402.4924000053



- To avoid an accidential tip ejection a clear haptic feedback is provided.
- The spring loaded tip cone improves user to user reproducibility and improved ergonomics by lowering tip attachment force
- The Reference 2 multi channel offers a high flexibility incl a spring-loaded tip cone which can be switched on/off optionally
- The secondary adjustment supports the easy adjustment of your Reference 2 for the most accurate pipetting of different liquids or other external conditions without the need for a full calibration.
- The quick lock enables an easy cleaning and maintenance of your pipette.
- The oustanding precision and accuracy of Reference 2 provides reliable results
- 4 digit volume display with magnifying window to allow for easier volume identification
- External edges made from stainless steel equip the pipette with outstanding robustness
- An embedded RFID chip contains all relevant data regarding the pipette and offers with Eppendorf TrackIT a simple identification and documentation