

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

Phone: 495 220 229, 495 220 394

Intramuscular Injection Training Buttocks

Order code: 4119.BTCSIM1



Information about product price on demand

Parameters

Cannula, injection, puncturing - filter

Buttock

ks

Quantitative unit

Quantitative unit

- The touch feeling of the model is similar to real human skin.
- The buttocks model is composed of epidermis layer, muscular layer and skeletal frame.
- The skin is made of soft silicon that leaves no needle marks.
- Ventrogluteal and dorsogluteal site injection practice is possible.
- Both side of buttocks can be used for practice.
- Drug injection and Automated drain
- Mechanical pump drains the injected liquid automatically.
- Self training and Evaluation mode are supplied.
- The location and depth of injection is detected and displayed.
- Animated video and Guidelines for injection procedures are supplied.
- The evaluation results are stored with Graphic database.
- Touch screen type user interface with tablet PC.
- Automated on-line connection with Bluetooth interface.
- User can select AC 110V/220V 60/50Hz power or Battery power.
- The efficacy of training can be maximized by graphic animations with anatomical representation.