



HELAGO-CZ, s.r.o.  
Commercial Register maintained by the Regional Court in Hradec  
Králové  
Section C, File 17879  
Kladská 1082  
500 03 Hradec Králové 3  
Company ID: 25 96 39 61, VAT: CZ 25963961  
Phone: 495 220 229, 495 220 394  
Fax: 495 220 154  
GSM gate: 602 123 096  
E-mail: [info@helago-cz.cz](mailto:info@helago-cz.cz)  
Web: <http://www.helago-cz.cz>

---

**B-510POL Trinocular Polarizing  
Microscope**

Order code: **0501.B510POL**



Cena bez DPH

3.167,00 Eur

Price with VAT

3.832,07 Eur

Parameters

Microscopes types

Laboratory

Serie

B-510

Number of eyepieces

Trinocular

Number of lenses

4

Zoom of lenses

4x, 10x, 20x, 40x

Illumination

LED

Quantitative unit

ks

**Observation Methods:** Brightfield, Transmitted polarization

**Head:** Trinocular, 30° inclined, rotatable 360°, ±5 dioptic adjustment on left eyepiece, interpupillary adjustment 50 – 75 mm, Fixed photo port, 50/50 ratio

**Eyepieces:** Plan wide field, PL 10x/22, high eyepoint. One of the eyepieces has a cross-hair reticle for optical centering

**Nosepiece:** Quintuple revolving nosepiece, rotation on ball bearings. Centering system for each objective.

**Objectives:** IOS W-PLAN POL 4x, 10x, 20x, 40x

**Stage:** Rotatable stage with locking mechanism and centering knobs. Vernier scale with accuracy 0.1 mm. Diameter 145 mm. Specimen slide clamps.

**Condenser:** Swing-out type condenser, sliding-in, N.A. 0.2/0.9, with centering system. With rotating polarizer filter.

**Bertrand Lens:**

- Swing-out type with centering mechanism for observation in conoscopy/orthoscopy. Rotatable analyzer from 0° to 90° with graduated scale. Tint plates included:
- 1° order red ( $\lambda$ )
- $\lambda/4$
- Quartz wedge

**Focusing:** Coaxial coarse and fine focusing mechanism (graduated, 0.002 mm) with upper stop. Adjustable tension of coarse focusing knob.

**Transmitted Illumination & Diaphragms:** High-efficiency X-LED3 (3.6 W), manual brightness control. Köhler illumination sys