

HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional Court in Hradec Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961
Phone: 495 220 229, 495 220 394

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz BT-CSCM - Male Catheterization Simulator

Order code: 4119.BTCSCM



Training environment more similar to clinical situation

Parameters

Quantitative unit

• We make a training environment which is more similar to clinical situation by using a wide range of catheters and a water-soluble lubricant jelly

ks

- Realistic resistance during insertion of urinary catheter
- When a catheter enters the bladder, artificial urine will come out
- External genitalia is separable from the buttocks/Easy to maintain
- Indwelling Catheterization and Intermittent Catheterization
- Practice with various size(Fr) of catheter

Features

- Life-like practical training
- Make position similar to the real patient
- Any water-soluble Jelly being used at clinical situation can be applied
- Training with various size of catheter is available
- Non-return valve : Prevent leaking after taking catheter out.
- Checking fluid after correct insertion of catheter at the right position.
- Foreskin (attachable / detachable)
- Genital organs and buttocks are attachable : easy to maintain

Specification

• Adult buttocks model: 3.5□ (430×280×200□)

• Storage case: 2.4[] (440×330×390[])

Total: Approximately 6□