



HELAGO-CZ, s.r.o.

Commercial Register maintained by the Regional Court in Hradec Králové

Section C, File 17879

Kladská 1082

500 03 Hradec Králové 3

Company ID: 25 96 39 61, VAT: CZ 25963961

Phone: 495 220 229, 495 220 394

Fax: 495 220 154

GSM gate: 602 123 096

E-mail: [info@helago-cz.cz](mailto:info@helago-cz.cz)

Web: <http://www.helago-cz.cz>

---

**Miele E 136 Insert**  
Order code: **1201.MIE69513601**



Information about product price on demand

Parameters

Quantitative unit

ks

- Suitable for 56 half petri dishes with Ø 100 mm
- equipped with 56 holders, height 70 mm, distance apart approx. 26 mm
- For use in lower basket

## Product details - E 136

Product category	
Washer-disinfectors, laboratory	•
Large-capacity washer-disinfectors, laboratory	•
Machine/accessory categorisation	
G 7825	•
PG 8527	•
PG 8535	•
PG 8536	•
PG 8583	•
PG 8583 CD	•
PG 8593	•
Product category	
Insert/module for petri dishes	•
Product properties	
Material	Stainless steel
Colour	Stainless steel
Distance between rails in mm	26
Rail height in mm	70
Area of application	
Reprocessing laboratory glassware	•
Capacity	
Petri dishes [number]	56
Dimensions and weight	
External dimensions, height in mm	485
External dimensions, width in mm	145
External dimensions, depth in mm	445
Net weight in kg	2.10