



HELAGO-CZ, s.r.o.  
Commercial Register maintained by the Regional Court in  
Hradec Králové  
Section C, File 17879  
Kladská 1082  
500 03 Hradec Králové 3  
Company ID: 25 96 39 61, VAT: CZ 25963961  
Phone: 495 220 229, 495 220 394  
Fax: 495 220 154  
GSM gate: 602 123 096  
E-mail: info@helago-cz.cz  
Web: http://www.helago-cz.cz

## **P93B-S - Catheterization Simulator Set BASIC**

Order code: **4111.1020842**



Cena bez DPH

708,00 Eur

Price with VAT

856,68 Eur

Parameters

With the P93 BASIC catheterization trainer set, both male and female bladder catheterization procedures can be realistically demonstrated, practiced and assessed. The easy-to-change genital inserted into the holder and held in place with magnets. Both inserts are made from soft material. Foreskin and labia are movable and the penis can be stretched. This way the student can practice all the necessary gestures for catheterization (e.g. disinfection). The material of the genital inserts gives a realistic feeling when inserting and removing the catheter into the bladder with realistic resistance and insertion depths.

**Direct learning:**

When the catheter has been correctly inserted, fluid runs out, just like with a real patient. The catheter can also be checked through the transparent bladder at any time.

**Various levels of difficulty:**

On the male insert, three different levels of narrowing of urethra can be set. In the narrowest setting, catheterization is not possible, leading to the conclusion that a suprapubic catheter is needed.

**All the benefits of the P93 BASIC catheterization simulators in the set at a glance:**

- Female and male catheterization with realistic resistance
- 3 levels of adjustable narrowing of urethra
- Soft and movable foreskin and labia
- Liquid outflow if catheterization is successfully carried out
- The situation can be checked through the transparent bladder
- Stands securely on the table thanks to the non-slip feet
- Magnetic connectors for quick set-up and dismantling
- Easy to clean, and the inserts can be completely removed

**Delivery contents:**

- Stands with magnetic connections for genital inserts
- Male and female catheterization set with transparent bladder
- Bottle and tube connectors

**Note:**

For catheterization, we recommend using catheters with only one or two drainage eyes opposite each other. Both silicone and PVC catheters can be used. We recommend the following catheters:

- Urethral narrowing "open" = 16 CH catheter or 14 CH balloon catheter
- Urethral narrowing "partial narrowing" = 14 CH catheter or 12 CH balloon catheter

**Catheters with a curved tip should not be used!****Tip:**

According to the guidelines for correct catheterization, even during the simulation, lubricating gel (not included) should be placed into the urethra to ensure that the catheter can be inserted smoothly. Alternatively, some lubricating gel can be placed on the tip of the catheter.