

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961

Phone: 495 220 229, 495 220 394

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz P94 - Birthing Stages Trainer Order code: **4111.1020628** 



Cena bez DPH

1.385,00 Eur

Price with VAT

1.675,85 Eur

Parameters

Obstetrics, gynecology - filter

Birthing simulators

Quantitative unit ks

• Detailed description

Trainees and students in midwifery and obstetrics can practice correctly assessing the progress of the birth by palpating, using the birth stages trainer. By practicing on the birth stages trainer, the standard of clinical practice in midwifery and obstetrics can be considerably improved.

The trainer consists of a set of 6 inserts that show the level of cervical dilation and effacing. The inserts can be used both as stand-alone trainers and as an extension module for birth simulator P90. The inserts look the same from the outside, but indicate, using a number code on the back, the condition of the cervix on the inside.

The following levels and birth stages are shown:

- Insert no. 1: Cervix closed (0.5 cm), level -2
- Insert no. 2: Cervix open 2 cm, level -2
- Insert no. 3: Cervix open 3 cm, half effaced, level -2
- Insert no. 4: Cervix open 5 cm, effaced, level -1
- Insert no. 5: Cervix open 7 cm, effaced, level -1
- Insert no. 6: Cervix open 9 cm, effaced, level 0

## **Delivery contents:**

- 6 birth stage inserts
- Safety plate with 4 screws (for the optional insert in birth simulator P90)
- Lubricant
- Carry bag

Do you want to print the page now? Hide some parts of page before printing



Print partially

Print the page as you see it now



Print complete page

(this box will not be printed)

Edit before printing

Hover your mouse over some part of page, a red rectangle will appear. Click it to hide the selected content of page.



Print!



Show hidden