



HELAGO-CZ, s.r.o.  
Commercial Register maintained by the Regional Court in  
Hradec Králové  
Section C, File 17879  
Kladská 1082  
500 03 Hradec Králové 3  
Company ID: 25 96 39 61, VAT: CZ 25963961  
Phone: 495 220 229, 495 220 394  
Fax: 495 220 154  
GSM gate: 602 123 096  
E-mail: info@helago-cz.cz  
Web: http://www.helago-cz.cz

**TruBaby X - Simulator of pediatric  
clinical skills**  
Order code: **4110.TB10001X**



Information about product price on demand

Parameters

Quantitative unit

ks

**The complete solution for pediatric clinical skills training**

**Airway Management**

- Lumbar Puncture

**PICC line insertion**

- Urethral Catheterization

**Needle Thoracentesis**

- Cardiopulmonary Resuscitation

**Peripheral Venous Cannulation**

- IO Tibia

*Why choose this model?*

- It's realistic! Weight, size and movement of a 50th Percentile 5 month old baby
- Ideal for practising DOP'S ( directly observed procedural skills)
- The model has been designed as fully interchangeable, all the parts are cost effective and can be replaced quickly when required.
- A Unique fully integrated fluid management system to replicate blood flow with no mess.
- Colour coded interface to allow easy prefill of fluids and air into the system
- The Arm and leg can detach to become separate task trainers for IV cannulation and IO procedures.
- Features the AirSim X airway which provides realistic and durable training and comes with a 5 year warranty.
- Delivered ready to use in durable carry case for easy transportation and safe storage

## Features and benefits

### Airway Management

- Realistic and durable Trucorp Airway
- Double naso-tracheal intubation
- Endotracheal tube insertion
- Bag-valve mask (BVM) ventilation techniques
- Full range of supraglottic device insertion
- Direct and video laryngoscopy
- Interchangeable head in case of repairs

### Peripheral Venous Cannulation (Arm)

- Contains the dorsal venous arch and metacarpal veins allowing needle cannulation at various locations
- Fully enclosed fluid management system giving realistic blood flashback and adjustable flow
- Fluids can be withdrawn and administrated
- 300+ needle incisions with no leaks using recommend needle gauge sizes
- The arm can become detachable and used as a stand alone hand IV task trainer for cannulation

### Peripheral Venous Cannulation (Foot)

- Contains the anterior arch vein, short and Great Saphonous veins and the superficial venous arch veins allowing needle cannulation at various locations
- Fully enclosed fluid management system giving realistic blood flashback and adjustable flow
- Fluids can be withdrawn and administrated
- 300+ needle incisions with no leaks using recommend needle gauge sizes
- The leg can become detachable and used as a stand alone Foot IV task trainer for cannulation

### Lumbur Puncture

- Ability to manoeuvre the baby into the required lateral decubitus position or sit upright
- Palpable vertebrae landmarks including the iliac crest
- Insertion sites include L3-L4 and L4-L5 spine locations
- Accurate needle placement allows for positive response and collection of simulated cerebrospinal fluid (CSF)
- Fluid is isolated and can be quickly refilled for a unique training experience
- Insert can last for 100+ incisions and is fully interchangeable

### PICC line insertion

- Complete PICC line procedure including guidewire insertion, dilation, and full cannulation
- Peripherally Inserted Catheterization either through the Basilic, Brachial, or Cephalic Veins, all visually located.
- Fully enclosed fluid management system giving realistic blood flashback and flow
- 300+ needle incisions with no leaks using recommend needle gauge sizes

### Needle Thoracentesis for a tension pneumothorax (2nd and 5th intercostal space)

- Realistic rib structures including xiphoid process and clavicle
- Palpable landmarks in the 2nd intercostal space midclavicular line and the 5th intercostal space mid auxiliary line
- Easily interchangeable inserts
- Each insert can last 150+ needle penetrations
- Realistic 'hiss' upon air escaping

**IO Tibia**

- Proximal Tibia interosseous needle insertion site
- Contains tibial tuberosity and patella anatomy
- Realistic feel when penetrating the medullary cavity
- Isolated insert with pre filled blood fluid for correct placement of IO needle
- Each insert will last 2-3 times due to destructive nature of the procedure
- The leg can become detachable and used as a stand alone IO task trainer in addition to the foot being used for IV cannulation

**Urethral Catheterization**

- Realistic anatomy to practice inserting catheter tube to extract urine
- Realistic Urethral and Vaginal opening
- Positive fluid flow when procedure has been successfully performed.
- Easy pre-fill mechanism

**CPR**

- Realistic rib structure with xiphoid process
- Full recommended depth of 1.5 inches can be achieved with full chest re coil

Weight: Approx. 5.5Kg

Dimensions: 62cm x 22cm x 12cm