




HELAGO-CZ, s.r.o.  
Commercial Register maintained by the Regional Court in Hradec  
Králové  
Section C, File 17879  
Kladská 1082  
500 03 Hradec Králové 3  
Company ID: 25 96 39 61, VAT: CZ 25963961  
Phone: 495 220 229, 495 220 394  
Fax: 495 220 154  
GSM gate: 602 123 096  
E-mail: info@helago-cz.cz  
Web: http://www.helago-cz.cz

**SURGICAL WOUND, WITHOUT  
SUTURE, 20CM**  
Order code: **4003.PW009A**



Cena bez DPH	235,00 Eur
Price with VAT	284,35 Eur

Parameters	
Quantitative unit	ks

This moulage is made of skin-friendly material and adheres by itself on the skin of a standardized patient, a manikin or on the supplied base. The moulage can be attached to any location of the body, for example over the knee to simulate a wound after a knee replacement or on the belly to simulate a cesarean section. The soft material moves with the body and can be attached over joints. If the moulage does not adhere well anymore, simply clean the backside with alcohol and it will stick strongly again. If the moulage shall be used for a longer course, we recommend a special skin glue.