



HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional Court in
Hradec Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961
Phone: 495 220 229, 495 220 394
Fax: 495 220 154
GSM gate: 602 123 096
E-mail: info@helago-cz.cz
Web: http://www.helago-cz.cz

BT-SEEM AIR - CPR & Airway Management Training Model

Order code: 4119.BTSEEMAIR



Cena bez DPH	2.662,00 Eur
Price with VAT	3.221,02 Eur

Parameters

Construction	Adult torso
Quantitative unit	ks

Features

- Training Chest Compression, Airway management, Ventilation for basic life support
- Intubation by using laryngoscope
- Insert the oral airway (OPA) and nasal airway(NPA) to the manikin
- By attaching offering magnetic stickers to AED trainer pads, AED trainer pads can be used semi-permanently
- Feedback about chest compression is being displayed through LED panel at the abdomen
- Training and evaluation of CPR with S/W
- Available to check evaluation result and print it out with S/W
- Self-training of CPR with a tablet or laptop / Group training with up to 30 manikins at the same time

Specification

Airway management

- Head tilt -Chin lift, Jaw thrust maneuver
- Ventilation with a bag-valve-mask
- Insert OPA/NPA
- Intubation with a laryngoscope
- Intubation with an endotracheal tube (ETT), laryngeal mask (LMA/l-gel), combi tube (Combi-Tube)
- Observable chest rising during the ventilation
- Click sound for the front teeth due to the excessive use of laryngoscope

Chest compression

* Wireless connection between manikin and tablet/smartphone via Bluetooth

- Easy-to-use and user-friendly S/W
- Possible to evaluate up to 30 manikins at the same time; tablet, printer and laptop are optional
- With its S/W, change the guide line before training
- With its S/W, display the compression depth, compression rate, hand position and also the ventilation volume
- With its S/W, check the result after training and save/print it
- Export the result data (Result in statistics) to excel format
- Score the training result and rank the saved scores
- Choose the mode between Training mode and Evaluation mode

*LED Display embedded in the manikin

- Real-time feedback of chest compression depth and rate with LED lamp
- Display total number of compression and the number of good compression
- Display the instant depth of compression
- Detect incomplete recoil
- Detect the good chest compression when depth and rate are good and give sound feedback (Sound On/Off)
- Display the feedback about ventilation such as volume and expiratory time
- Display the total number of ventilation and the number of good ventilation
- Display the ventilation volume

* Manikin body

- Observable chest rising during ventilation
- By attaching offering magnetic stickers to AED trainer pads, AED trainer pads can be used semi-permanently
- The storage bag can be used as a training mattress
- Use batteries / AC adaptor

Components

- 1 Adult Torso manikin
- 1 Storage case
- 1 Lubricant for intubation
- 1 AC Power Adaptor, User Manual