



HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional Court in
Hradec Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961
Phone: 495 220 229, 495 220 394
Fax: 495 220 154
GSM gate: 602 123 096
E-mail: info@helago-cz.cz
Web: http://www.helago-cz.cz

DW-HL340 Ultra Low Freezer For Laboratory And Medical Freezer

Order code: **2203.DWHL340**



Information about product price on demand

Parameters

Device type	Freezer
Description_	4 wheels
Volume brutto	340 l
The temperature range in the freezer	-40 až -86°C
Number of sockets	3 shelves (4 inner doors)
Defrosting	manual
Refrigerant	mixed gas
Quantitative unit	ks

The digital temperature display can show various parameters simultaneously, such as temperature inside the cabinet, power voltage and environmental temperature, indicating operating status clearly The high-precision microcomputer temperature control system The well-developed alarm system has multiple alarm functions: high temperature/low temperature alarm, door ajar alarm, power failure alarm, sensor failure alarm, low battery alarm, alarm of poor heat dissipation of condenser, system failure alarm, etc. The keyboard lock function and password protection function can prevent any adjustment of operating parameters without permission The rotating handle facilitates door opening The adjustable layered structure is convenient for storage of various things