




HELAGO-CZ, s.r.o.  
Commercial Register maintained by the Regional  
Court in Hradec Králové  
Section C, File 17879  
Kladská 1082  
500 03 Hradec Králové 3  
Company ID: 25 96 39 61, VAT: CZ 25963961  
Phone: 495 220 229, 495 220 394  
Fax: 495 220 154  
GSM gate: 602 123 096  
E-mail: info@helago-cz.cz  
Web: <http://www.helago-cz.cz>

## **LPBP-40-L - LumbarPunctureBaby System Training Package - Light Tone**

Order code: **4129.LPBP40L**



Information about product price on demand

Parameters

Cannula, injection, puncturing - filter

Puncture, anesthesia, injections

Simulab's LumbarPunctureBaby System Training Package simulates a two-week-old infant that can be positioned either in the lateral decubitus position or held upright in the sitting position. The body form is anatomically correct with an iliac crest and umbilicus. The replaceable tissue has L3 - L5 vertebrae with a partial and gluteal fold.

**Includes:**

- LumbarPunctureBaby Body Form
- Two Replaceable Tissues
- Molded Body Form Tray
- Carrying Case
- Simulated IV Bag
- Table Top IV Pole
- Bottle of Red Concentrate

**Features:**

- Anatomically correct two-week old infant with flexible body form
- Ultrasound compatible
- Visible and palpable landmarks include umbilicus, gluteal fold, iliac crest and vertebrae
- Insertion sites include L3-L4 and L4-L5
- The flexible body adds realism when flexing the infant, simulating moving the interspinous process from a neutral to open position
- Accurate needle placement allows for positive response and collection of simulated cerebrospinal fluid (CSF)
- Simulated epidural venous plexus and bony spinous process provides the user with feedback for improper needle placement
- IV Bag provides control of the drip rate at the procedure site

**Skills:**

- Use of the lateral decubitus and sitting positions for pediatric lumbar puncture
- Palpate external landmarks
- Proper technique for lumbar puncture (spinal tap) procedure