



HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional Court
in Hradec Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961
Phone: 495 220 229, 495 220 394
Fax: 495 220 154
GSM gate: 602 123 096
E-mail: info@helago-cz.cz
Web: <http://www.helago-cz.cz>

**LP-30-L - Lumbar Puncture Trainer with
Simulated IV Bag (IVB-20) - Light Tone**

Order code: **4129.LP30L**



Information about product price on demand

Parameters

Cannula, injection, puncturing - filter

Quantitative unit

Puncture, anesthesia, injections

ks

Safely master critical techniques for lumbar puncture insertion with Simulab's easy-to-use, ultrasound-guided **Lumbar Puncture Trainer**.

Our true-to-life simulator sits upright or in a lateral decubitus position and features the lumbar vertebrae, iliac crest, spinous process, ligamentum flavum, and epidural space; with a fluid-filled compartment for simulated cerebrospinal fluid. The lumbar puncture allows the user to penetrate the subarachnoid space and CSF is present when the stylet is removed.

The trainer supports three varieties of lifelike, ultrasound-compatible replaceable tissues that mirror significant anatomical differences found in obese, geriatric, and normal patient populations. Every tissue allows users to experience variances in disc thickness and soft tissue depth on top of the vertebrae.

Included

- Simulated IV Bag
- Normal Tissue

Features

- Self-healing replaceable tissue is highly durable for repeated procedures
- Replaceable tissue variations allow for progressive skill development of diverse patient populations
- Durable and replaceable tissue
- Most realistic tactile feedback on the market

Skills

- Use ultrasound or palpation to identify key landmarks
- Experience lifelike needle resistance, including pops when needle transverses ligamentum flavum and dura
- Allows students to collect CSF
- Practice the procedure in the left lateral decubitus or sitting position