




HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional Court in
Hradec Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961
Phone: 495 220 229, 495 220 394
Fax: 495 220 154
GSM gate: 602 123 096
E-mail: info@helago-cz.cz
Web: <http://www.helago-cz.cz>

PAC-20-L - Ultrasound Paracentesis
Model - Light Tone
Order code: **4129.PAC20L**



Information about product price on demand

Parameters

Phantoms and simulators

Ultrasound phantoms

Simulab's Ultrasound Paracentesis Trainer ultrasound compatible trainer allows the user to perform diagnostic and/or therapeutic paracentesis. The anatomically correct trainer is positioned to simulate a patient sitting up at a 45-degree angle with a distended abdomen and includes the pubic symphysis, iliac crest, and umbilicus. The procedure can be performed at the midline below the umbilicus or 2 fingerbreadths anterior and 2 fingerbreadths cephalad to the anterosuperior iliac spine. The inferior epigastric vessels can be visualized under ultrasound.

Includes:

- Paracentesis Base
- Replaceable Tissue (PACT-20)
- 7/16" Nut Driver

Features:

- Ultrasound compatible with replaceable tissue
- Internal anatomy includes superficial epigastric vessels, liver and spleen, rectus abdominal muscles, and intestines
- Allows up to one liter of intraperitoneal fluid removal
- Replaceable tissue is durable and allows for repeated use
- Palpable anatomy and realistic needle response

Skills:

Use an ultrasound probe to minimize risks and learn to:

- Locate the pocket of ascites
- Gauge the distance from skin to fluid, and avoid bowel
- Look for blood vessels beneath the skin using a vascular probe
- Drain the fluid from the peritoneal cavity
- Use the traditional lateral gutter technique or midline linea alba technique
- Palpate anatomic landmarks significant to the procedure