

**P RTP-30-L - Perineal Repair Training
Package - Light Tone**
Order code: **4129.PRTP30L**



Information about product price on demand

Parameters

Obstetrics, gynecology - filter

Common gynecological and obstetric actions

Simulab's Perineal Repair Training Package is an affordable, realistic, and easy to use solution for episiotomy and perineal repair. The package includes the training base and two replaceable tissues.

The soft-tissue model provides a highly realistic and practical solution for training healthcare professionals in repairing first to 4th-degree perineal tears, labial lacerations, catheter insertion, episiotomy, and end-to-end or overlap repair.

The trainer includes the vaginal, anal canals, urethra, rectovaginal fascia, external sphincter, internal sphincter, skin and muscle layers, and vaginal and rectal mucosa. The soft-tissue model allows for incisions or simulated tears along all sides of the vaginal canal, including midline, mediolateral, and lateral.

The replaceable tissue allows for multiple cuts and repeatable realistic suturing of the mucosa, muscles, skin, and subcutaneous layers. The urethral canal allows for catheter insertion to repair simulated tears along the top of the vaginal canal near the urethral canal.

The tray form overlays the tissue and provides tension that enhances the simulation. A bony anatomical landmark on the frame simulates the ischial spine allowing the user to perform a pudendal block. An adjustable distention rod adds tension to provide more realistic access to deeper tissues during 3rd and 4th-degree repairs.

Anatomy Includes:

- Vaginal & Anal Canals
- Urethra
- Labia
- Vaginal Mucosa
- Rectovaginal Fascia
- Transverse Perineal Muscle
- External & Internal Sphincter
- Ischial Spine

Includes:

- Perineal Repair Trainer Base
- Two Perineal Repair Replacement Tissue

Features:

- Replaceable tissue allows for multiple incision sites and suturing
- Realistic tissue properties for training 4th-degree repair that includes a retractable sphincter and rectal mucosa
- Adjustable tension for variety in training situations

Skills:

- Incision training for midline, mediolateral and lateral incisions
- Insertion of a vaginal retractor and retraction with the non-dominant hand
- Suturing of realistic tissue layers
- Labial repair
- Perform a 1st to 4th-degree repair
- Practice urethra catheter insertion