

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3

Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz **BDT-10 - Bean Drop Task**

Order code: **4129.BDT10**







Information about product price on demand

Training 2 consecutive Teeting 2 of 3 consecutive

24

SIMULAB

Parameters

Quantitative unit ks

A proficiency-based task trainer, designed to teach the technical skills required in basic laparoscopic surgery. The Bean Drop Task requires the user to grasp and move each bean 15 cm from the Saucer to the Final Elevated Location Cup. Use with any of Simulab's Laparoscopic Trainers, or secure to table for classroom, clinic, home, or office use.

Features:

• This task includes a Right and a Left-Handed Template along with proficiency metrics for self-assessment or skills lab.