




HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional Court in
Hradec Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961
Phone: 495 220 229, 495 220 394
Fax: 495 220 154
GSM gate: 602 123 096
E-mail: info@helago-cz.cz
Web: http://www.helago-cz.cz

LF01121 - Advanced Venipuncture and Injection Arm

Order code: 4107.LF01121



Cena bez DPH

1.066,00 Eur

Price with VAT

1.289,86 Eur

Parameters

Cannula, injection, puncturing - filter

Horní končetiny

Quantitative unit

ks

This revolutionary training arm provides complete venous access for IV therapy and phlebotomy, plus sites for intramuscular and intradermal injections. An extensive 8-line vascular system allows students to practice venipuncture at all primary and secondary locations, including starting IVs and introducing Over the Needle IV Catheters. The venous system simplifies setup with only one external fluid bag supplying artificial blood to all veins simultaneously. The dorsal surface of the hand includes injectable metacarpal, digital, and thumb veins. The antecubital fossa includes the median cephalic, median basilic, and median cubital veins. Venipuncture can also be performed along the basilic, cephalic, accessory cephalic, and median antebrachial veins. Intramuscular injections may be performed in the deltoid muscle and intradermal injection sites located in the upper arm. Intramuscular injections into the deltoid muscle are enhanced by the soft, lifelike skin and by the natural bony landmarks in the region. Intradermal injections using distilled water will create characteristic skin welts at designated sites on the upper arm.

The realism of the Life/form® Advanced Venipuncture and Injection Arm is truly amazing. The soft, flexible fingers are molded separately with extreme attention paid to every detail — right down to the fingerprints! Flexion of the wrist helps students develop manipulation skills. The replaceable skin rolls as the veins are palpated, and a discernable “pop” is felt when entering the veins. The Life/form® molding process reproduces the fine details of the skin to make the arm look and feel alive. Valves in the veins can be seen and palpated on the skin surface.

Simulated veins and skins are completely replaceable to keep this training arm looking and working like new. Under normal use, hundreds of injections may be performed before the veins or skins need to be replaced. Complete replacement kits are available and easy to use. To

extend the life of the veins, an aerosol sealant is available to seal punctures and prevent leakage.

Includes 2 qts. simulated blood, 3 cc syringe, 12 cc syringe, needle, two fluid bags, instruction manual, and carrying case. Three-year warranty. Sh. wt. 22 lbs. (9.98 kg).

- LF01126 - Advanced Venipuncture and Injection Arm - Black

Available Supplies

- LF01122 - Skin and Vein Replacement Kits - White
- LF01123 - Skin and Vein Replacement Kits - Black
- LF01124 - Vein Replacement Kit
- LF01022 - Fluid Supply Stand
- LF01130 - Fluid Supply Bag
- LF00845 - Blood - Quart - Sh. wt. 0.38 lbs. (0.17 kg)
- LF00846 - Blood - Gallon - Sh. wt. 0.88 lbs. (0.40 kg)
- W09919 - REN Cleaner