


**50124 - Vermiform Appendix -  
 Retrocaecal**  
 Order code: **4104.50124**



Cena bez DPH	142,00 Eur
Price with VAT	171,82 Eur

Parameters	
Surgery - filter	Basic surgery
Quantitative unit	ks

For practising the laparoscopic removal of the vermiform appendix. There are 3 variations of this model which present the appendix in the following positions:

- normal (as found in 32% of patients)
- post-ileal (as found in 1% of patients)
- retrocaecal (as found in 64% of patients)

Each model has the additional feature of a simulated peptic ulcer, for closure using peritoneum and mesentery.

## Skills

- Excision and division of peritoneum
- Identification and exposure of appendix
- Mobilisation and division of vessels
- Removal of appendix
- Inspection of stump of appendix
- Repair and closure of peptic ulcer

## Features

- Highly realistic model of appendix, caecum, and ileum stump presented anatomically
- The variations available in the range of appendices provide increased levels of difficulty
- Convincing fluid filled vessels, caecum contents, and tissues encourage correct handling of tissue
- Vessel incorporates a luer lock to which the Mock Blood Giving Set may be attached to provide fluid flow for the appendicular artery.  
Use Mock Blood -Arterial (250 ml) - **60653**
- Ideal for training course requirements; the model is designed for swift attachment to and removal from the Soft Tissue Retaining Set