

PM 200 - Planetary Ball Mill

Order code: **1306.206400001**



Information about product price on demand

Parameters

Description

Převodový poměr 1:-2

Quantitative unit

ks

- extreme speeds for particularly high fineness down to submicron range
- 2 grinding stations
- grinding jars with volumes from 12 ml to 125 ml, in 6 materials
- suitable for long-term trials and continuous use
- reproducible results due to energy and speed control
- direction reversal
- programmable starting time
- power failure backup ensures storage of remaining grinding time
- 1-button operation with graphic display
- 10 combinations of grinding parameters can be stored
- measurement of input energy helps to optimize grinding parameters
- automatic grinding chamber ventilation for cooling the grinding jars
- 2-year warranty, CE-conforming