

PP00031 - Traction Splint Trainer

Order code: **4107.PP00031**



Cena bez DPH

3.455,00 Eur

Price with VAT

4.180,55 Eur

Parameters

Quantitative unit

ks

Traction splint application requires the ability to regulate how much force is applied to a patient's leg. There has never been a way to monitor that, until now. The Traction Splint Trainer can be set to correspond to an individual's body weight. This allows the instructor to determine if the student is correctly applying the rule of thumb that states "10% of the body weight traction results in adequate relief for, but no further damage, to the patient." Adjust for body weights ranging from 45 to 135 kg. The electrical, push-button panel makes simplified selection easy. This AC electrical trainer comes completely assembled and requires attachment to a base to secure it during training sessions. After applying traction, the deformity of the thigh recedes and splint can be applied, including any standard ankle hitch.