



Information about product price on demand

Parameters

Quantitative unit

ks

Model has 14 or 18 parts, can be supplied with transparent plastic falx cerebri and tentorium cerebelli. The location of cavities cerebrospinal dura mater is marked on falx cerebri. Pacchioni version is shown in the upper sagittal sinuses. Track venous blood and cerebrospinal fluid is easily recognizable. Anterior falx cerebri is attached to the crista galli and falx cerebri at the upper parietal bones. The model is easily severable and shows the position of the sinus dura scleriosis (blue).

Weight: 66 g