



HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional Court in Hradec Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961
Phone: 495 220 229, 495 220 394
Fax: 495 220 154
GSM gate: 602 123 096
E-mail: info@helago-cz.cz
Web: http://www.helago-cz.cz

Sakamoto Suction Simulator

Order code: **4109.M175**



Cena bez DPH

2.276,00 Eur

Price with VAT

2.753,96 Eur

Parameters

Quantitative unit

ks

Basis of respiratory management Endotracheal intubation can be learned.

This model is a simulator with a lifelike texture, and various suction methods can be learned.

The suction of secretions not only in the oral/nasal cavities but also in tubes inserted via the oral cavities is possible.

Since the trachea has already been incised, suction with tracheal cannula insertion can also be performed.

Features

The head can be moved to the left and right.

Practice

Various types of suction training are possible

Suction of secretions in the oral cavity

Suction of secretions in inserted tubes

Suction of secretions in the nasal cavity

Suction of secretions in the tracheal cannula

Washing of the bronchi

Washing of the oral cavity