




HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional Court in Hradec
Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961
Phone: 495 220 229, 495 220 394
Fax: 495 220 154
GSM gate: 602 123 096
E-mail: info@helago-cz.cz
Web: http://www.helago-cz.cz

W19348 - Pressurised Cannulation Pad

Order code: **4111.1005167**



Information about product price on demand

Parameters

Cannula, injection, puncturing - filter

Kanyly

Quantitative unit

ks

The Pressurised Cannulation Pad has been introduced to provide a portable intravenous trainer and act as an alternative to using the traditional arm. These trainers are realistic and achieve relatively low running costs.

- There is one blood vessel that is set into a fleshy pad at variable depths
- The Pressure is provided by a simple squeeze on the bulb into a fluid reservoir
- Both skin pad and vessels are replaceable

Measurements: 27 × 24,5 × 10 cm

Weight: 1,0 kg

Options and Replacements:

- **W19349 - Cannulation Pad**
Measurements: 27 × 24,5 × 10 cm
Weight: 800 g
- **W19350 - Replacement Tubing**
Measurements: 20 × 1,2 × 1 cm
Weight: 100 g