



HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional Court in Hradec
Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961
Phone: 495 220 229, 495 220 394
Fax: 495 220 154
GSM gate: 602 123 096
E-mail: info@helago-cz.cz
Web: http://www.helago-cz.cz

SB46506 - "Pat" Pressure Ulcer Model

Order code: **4107.SB46506**



Cena bez DPH

907,00 Eur

Price with VAT

1.097,47 Eur

Parameters

Quantitative unit

ks

A unique, compact, comprehensive, and realistic model that is an effective tool for identifying and staging pressure ulcers. Also a great visual aid to show what can occur without proper care. Displays Stage I, II, III, and IV pressure ulcers, an unstageable eschar/slough wound, an intact blister, open wound, undermining, tunneling, and exposed bone, tendon, and muscle. The darkly pigmented skin section shows

how different a Stage I, blood blister, and suspected Deep Tissue Injury may appear depending on the pigmentation. The DTIs have a “spongy/soft” feel, which is especially helpful in darkly pigmented patients where visualization of the injury may not be apparent. Routine cleansing and dressing changes can be taught and practiced on all wounds, and you won’t need to constantly cleanse the skin of residual adhesive.