

**LF03674 - Child CRiSis™ Starter
Torso**

Order code: **4107.LF03674**



Cena bez DPH

3.450,00 Eur

Price with VAT

4.174,50 Eur

Parameters

Construction

Adult torso

Quantitative unit

ks

- Airway Management
- Arrhythmia Recognition

- Defibrillation
- ECG Monitoring (4 sites)
- Manual Carotid Pulse
- Pacing

The Life/form® Child CRiSis™ Starter Torso can be completely upgraded to a dramatic, full PALS training manikin. Present your students with the challenges they may face in the real world. This airway management trainer is perfect for practicing skills on pediatric patients and represents the size of an 8-year-old patient. Practice intubation, ventilation, suction, CPR, and jaw thrust techniques. Manual carotid pulse. Realistic anatomy and landmarks include teeth, tongue, oral and nasal pharynx, larynx, epiglottis, arytenoids, false and true vocal cords, trachea, lungs, esophagus, and stomach. The trainer allows you to practice oral, digital, and nasal intubation, as well as E.T., E.O.A., P.T.L., L.M.A., E.G.T.A., Combitube®, and KING System insertion. Separate lungs for auscultation. Inflatable stomach bladder indicates esophageal insufflation. With its slightly anterior position, the Life/form® Child CRiSis™ Starter Torso is ideal for ACLS, nursing, paramedic, and EMT training at every level. This trainer with defibrillation chest skin lets you practice defibrillation using standard manual, automatic, or semi-automatic external monitor defibrillators. An internal load box absorbs the full strength of every shock to protect students and equipment. It is also possible to monitor the manikin, just like a real patient, and it is compatible with all standard brands and types of defibrillators, monitors, and patient simulators. Features four ECG sites and two defibrillation sites. This torso is designed to be the foundation of the Life/form® Child CRiSis™ Manikin. Arms and legs can be added later to provide IV access, intraosseous infusions, and blood pressure monitoring. If you do not have an ECG simulator and want to practice ECG recognition, you may want to consider adding the Life/form® Interactive ECG Simulator to your training (LF03670). If advanced airway management is preferred, consider starting with the Life/form® Advanced Child CPR/Airway Management Torso with Defibrillation Features (LF03763). Pump spray lubricant and instruction manual are included.