

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961

Phone: 495 220 229, 495 220 394 Fax: 495 220 154

GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

## Eppendorf Reference® 2, variable 0.5 - 5 ml

Order code: 1402.4924000100





Cena bez DPH
Price with VAT

**Parameters** 

Pipettes S nastavitelným objemem

Description pro špičky 0,1 - 5ml, 120 mm

Quantitative unit ks

• The unique one-button operation offers a fast and ergonomic handling with reduced operating effort and active aerosol reduction.

To avoid an accidential tip ejection a clear haptic feedback is provided.

462,22 Eur

- The spring loaded tip cone improves user to user reproducibility and improved ergonomics by lowering tip attachment force
- The Reference 2 multi channel offers a high flexibility incl a spring-loaded tip cone which can be switched on/off optionally
- The secondary adjustment supports the easy adjustment of your Reference 2 for the most accurate pipetting of different liquids or other external conditions without the need for a full calibration.
- The quick lock enables an easy cleaning and maintenance of your pipette.
- The oustanding precision and accuracy of Reference 2 provides reliable results
- 4 digit volume display with magnifying window to allow for easier volume identification
- External edges made from stainless steel equip the pipette with outstanding robustness

•	An embedded RFID and documentation	chip contains all releva	nt data regarding the	e pipette and offers	with Eppendorf TrackIT	a simple identification