



HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional Court in
Hradec Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961
Phone: 495 220 229, 495 220 394
Fax: 495 220 154
GSM gate: 602 123 096
E-mail: info@helago-cz.cz
Web: http://www.helago-cz.cz

US-16 - Intravesical Urine Volume Measurement Training Simulator

Order code: **4103.US16**



Information about product price on demand

Parameters

Quantitative unit

ks

Pocket-sized ultrasound is becoming popular and start to be utilized in adult day care to check intravesical urine volume instead of catheterization which may cause urinary infection.

Training skills

- Handling and manipulation of transducers on bladder ultrasound
- Scanning the bladder to measure urine volume and findings
 - 4 variations of interchangeable inserts feature different patient scenarios
 - Urine volume 50 ml
 - Urine volume 150 ml
 - Urine volume 300 ml
 - Urinary retention + balloon catheter

Set Includes:

- 1 pelvic model
- 1 carrying container

Specifications:

- Body size: approx. W30×D26×H18cm
- Phantom size: approx. W14×D12×H13cm
- Body weight: approx. 3.2kg
- Phantom weight: approx.1.5kg