



HELAGO-CZ, s.r.o.

Commercial Register maintained by the Regional Court in Hradec Králové

Section C, File 17879

Kladská 1082

500 03 Hradec Králové 3

Company ID: 25 96 39 61, VAT: CZ 25963961

Phone: 495 220 229, 495 220 394

Fax: 495 220 154

GSM gate: 602 123 096

E-mail: [info@helago-cz.cz](mailto:info@helago-cz.cz)

Web: <http://www.helago-cz.cz>

---

**Miele O 175 Upper basket**

Order code: **1201.MIE69117501**



Information about product price on demand

Parameters

Quantitative unit

ks

- Suitable for narrow-necked and laboratory glassware
- Vertical clearance: 170 mm<sup>1</sup>
- Equipped with connection for hot air drying
- Powder dispensing not possible

## Product details - O 175

<b>Product category</b>	
Washer-disinfectors, laboratory	•
<b>Machine/accessory categorisation</b>	
PG 8535	•
PG 8536	•
<b>Product category</b>	
Insert/module for narrow-necked glassware	•
Upper basket	•
<b>Product properties</b>	
Material	Stainless steel
Colour	Stainless steel
Maximum vertical clearance in mm	170
<b>Area of application</b>	
Reprocessing laboratory glassware	•
<b>Capacity</b>	
Narrow-necked glassware [number]	33
<b>Dimensions and weight</b>	
External dimensions, height in mm	412
External dimensions, width in mm	475
External dimensions, depth in mm	531
Net weight in kg	6.30
<b>Standard accessories</b>	
E 351 injector nozzles Ø 4 mm, length 160 mm [number]	33
E 353 clips for injector nozzles E 351 [number]	33