

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394 Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

Rotor A-4-81, without buckets Order code: 1309.5810743001



Cena bez DPH Price with VAT

Parameters Device type 1.494,00 Eur

1.807,74 Eur

Accessories Quantitative unit ks

- Swing-bucket rotor with max. speed 3,220 × g (4,000 rpm)
- Rotor A-4-81 accommodates tubes and bottles from 0.2 mL to 500 mL as well as plates
- Aerosol-tight caps¹⁾ for safe centrifugation of hazardous samples optionally available
- Mixed loading with 2 plate buckets and 2 rectangular tube buckets allows for increased versatility
- Buckets, caps and adapters are autoclavable (121 °C, 20 min)
- ¹⁾ Aerosol-tightness tested by Public Health England, Porton Down, UK.