



HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional
Court in Hradec Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961
Phone: 495 220 229, 495 220 394
Fax: 495 220 154
GSM gate: 602 123 096
E-mail: info@helago-cz.cz
Web: http://www.helago-cz.cz

Rotor FA-6x250, for 6 x 250 mL bottles, incl. aerosol-tight Eppendorf QuickLock® rotor lid

Order code: **1309.5895175007**



Cena bez DPH

4.750,00 Eur

Price with VAT

5.747,50 Eur

Parameters

Device type

Accessories

Quantitative unit

ks

- Fixed-angle rotor for up to 6 x 250 mL bottles¹⁾
- Ideal for the pelleting of microorganisms and mammalian cells
- 12 different adapters for broad range of additional applications
- Max. speed 15,054 x g (10,100 rpm)
- Eppendorf QuickLock rotor lid for fast lid locking
- Aerosol-tight²⁾ for safe centrifugation of hazardous samples
- Low weight of only 5,450 g for stress-free work
- Anodized aluminum for high chemical resistance
- Rotor, lid and adapters are autoclavable to reduce risk of contamination
- See tab 'Additional information' for overview of adapters available
- ¹⁾ To run this rotor, a software update is required for earlier centrifuges. Please contact your local Eppendorf organization or dealer for further information. ²⁾ Aerosol-tightness tested by Public Health England, Porton Down, UK. Maximum 80% of maximum filling volume.