

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394 Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

## Rotor S-4x500, for 500 mL rectangular buckets or MTP/Flex buckets, incl. 4 x 500 mL rectangular buckets Order code: 1309.5895170005



Accessories

ks

Device type Quantitative unit

• Medium capacity swing-bucket rotor with max. speed 3,220 × g (4,000 rpm)

- Rotor S-4x500 accommodates tubes and bottles from 0.2 mL to 500 mL as well as plates
- Aerosol-tight caps for safe centrifugation of hazardous samples optionally available<sup>1)</sup>
- Mixed loading with 2 plate buckets and 2 rectangular tube buckets allows for increased versatility
- Buckets, caps and adapters are autoclavable
- <sup>1)</sup>Aerosol-tightness tested by Public Health England, Porton Down, UK.