

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961 Phone: 495 220 229, 495 220 394 Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz

TruIV Block Order code: 4110.TIV100



Cena bez DPH	697,00 Eur	Eur
Price with VAT	843,37 Eur	
Parameters		
Cannula, injection, puncturing - filter	Ca	annulas
Ouantitative unit	ks	5

Part of the TruUltra product range, the TruIV Block allows learners to gain proficiency in the skills associated with IV cannulation. Featuring 8 veins (superficial and deep) ranging from 4-8mm in diameter, users can practise the use of ultrasound during IV procedures.

Model Features

- Realistic features: Features 8 veins (4 superficial and 4 deep) plus fascia layers for enhanced realism.
- Target various veins ranging from 4mm to 8mm in diameter. Ultrasound imaging is incredibly lifelike
- · Confirm successful IV cannulation procedures by withdrawing blood
- Ready to begin training in under 5 minutes
- Unique, self-regenerating TruUltra materials diminish the appearance of needle tracks to facilitate a high volume of repeated practice. The insert will facilitate approx. 4000 needle penetrations before a replacement insert is required
- · Replacement insert is cost-effective and easy to replace with no disruption to teaching time
- Delivered in a compact, protective carrier case perfect for training in different locations
- Longitudinal and transverse anatomical viewing options
- Compatible with all ultrasound machine brands

Medical Procedure Training

- Ultrasound-guided IV cannulation
- Probe positioning & movement
- Recognition of veins in soft responsive tissue