

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3 Company ID: 25 96 39 61, VAT: CZ 25963961

Phone: 495 220 229, 495 220 394

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz **C0008 - Control the bleed**

Order code: **4130.C0008**



Cena bez DPH 1.036,00 Eur
Price with VAT 1.253,56 Eur

Parameters

Quantitative unit ks

Control the bleed task trainer. Realistic and functional. With adjustable difficulty.

Made out of prosthetic grade silicone and 3D printed parts. Can be used for both wet and dry training. The blood flow speed can be easily adjusted.

LEARNING OBJECTIVES

Use it to teach:

- wound packing
- hemostatic devices application
- tourniquet application
- improvised tourniquets
- itclamp application
- using a foley catheter to stop the bleeding in the neck area

Main Features

hand-made with attention to details

Each wound was carefully copied from real-life examples.

bleeding neck

compatible with ITClamps and Foley catheters

adjustable blood flow

The blood flow speed can be easily adjusted to the training needs.

stable base

We added a stable base that mimics a shoulder joint to make the training more realistic.

hand-made with attention to details

There are three bleeding points in the tourniquet application arm in the place of brachialis, superior ulnar collateral and profunda brachii arteries.

Box contents

- wound packing arm
- tourniquet arm
- bleeding neck
- fluid container & tubing
- absorbent cloth
- silicone pad
- USBwith instructions
- carry bag