

**ARTP-20 - Arthrocentesis Training
Package (ARC-20, ARCT-20)**

Order code: **4129.A RTP20**



Information about product price on demand

Parameters

Quantitative unit

ks

Practice diagnosing and managing the presence of knee effusion with the anatomically correct simulator in this Arthrocentesis Training Package. The trainer represents an extended left leg with ultrasound-compatible areas, including the patella, patellar ligament, tibia, fibula, femur, synovial sac, and synovial fluid. Synovial fluid can be aspirated from a joint cavity using the medial or lateral approach. Insertion sites include suprapatellar and parapatellar. For added package savings, the Arthrocentesis Training Package comes with two replaceable tissues.

Includes:

- Arthrocentesis Knee Model
- Two Replaceable Arthrocentesis Tissues ARCT-20 One 60 ml Luer-lock syringe

Features:

- Ultrasound compatible
- Highly durable replaceable tissue for multiple uses
- Realistic tactile feedback includes:
 - A sensation of bony contact if the needle hits the patella
 - A sensation of bony contact if the needle hits the femur
- Inability to aspirate the syringe while the needle tip is in the soft tissue superficial to the joint capsule
- Easy aspiration of joint fluid into the syringe once entry is achieved
- Fluid can be left clear or be colored if desired
- Ability to increase or decrease the size of effusion with up to 60 ccs of fluid

Skills:

- Learn to ballot or milk the suprapatellar pouch
- Needle aspiration and injection techniques
- Palpating anatomic landmarks significant to the procedure
- Realistic tactile feedback throughout the procedure