




HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional Court in Hradec
Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961
Phone: 495 220 229, 495 220 394
Fax: 495 220 154
GSM gate: 602 123 096
E-mail: info@helago-cz.cz
Web: http://www.helago-cz.cz

SNMA-21 - SonoMan AAA Training Module v2.1

Order code: 4129.SNMA21



Information about product price on demand

Parameters

Phantoms and simulators

Ultrasound phantoms

The AAA Module for the SonoMan® System Diagnostic Ultrasound Trainer is an extension to the SonoMan System FAST Module.

The AAA Module allows students to evaluate and confirm a high-risk diagnosis with a patient presenting with a life-threatening abdominal aortic aneurysm. The AAA Module includes a mentor, student, and instructor mode. In Mentor mode, the user has access to an annotated version of each image. Instructor mode allows the instructor to include additional annotations for advanced learning. The image and video library include representations of two normal patients and three abnormal patients where abnormalities can be diagnosed.

Note:

- This is a software module for the SonoMan Diagnostic Ultrasound Trainer (sold separately) Contact your Account Representative to learn more.
- PC not included

Includes:

- AAA Software Module for SonoMan System

Features:

- Landmarks on the body used for identification of window locations for the exam
- Seven targets in the window
- Transverse and longitudinal viewing
- All images and videos are of actual patients - both normal and abnormal M
- Mentor mode allows users to activate the onscreen labels of any of the digital images
- Each module includes an instructor (teaching) and student (exam) mode

Skills:

- Read diagnostic ultrasound imaging
- Reading and understanding the various normal and abnormal views
- Identify external and internal landmarks for a AAA exam