

HELAGO-CZ, s.r.o. Commercial Register maintained by the Regional Court in Hradec Králové Section C, File 17879 Kladská 1082 500 03 Hradec Králové 3

Company ID: 25 96 39 61, VAT: CZ 25963961

Phone: 495 220 229, 495 220 394

Fax: 495 220 154 GSM gate: 602 123 096 E-mail: info@helago-cz.cz Web: http://www.helago-cz.cz 1024776 - Practi-3% Hypertonic Saline 250 mL I.V. Solution Bag (×1)

Order code: 4132.1024776



Cena bez DPH

3,60 Eur

Price with VAT

4,36 Eur

Parameters

Quantitative unit ks

The Practi-3% Hypertonic Saline 250 mL I.V. Solution Bag for clinical training, simulates a 3% hypertonic saline 250 mL I.V. Solution Bag, often used in Tactical Combat Casualty Care (TCCC.). The bag contains distilled water, making it suitable for use on manikins/skill trainers that are compatible with fluids.

The IV-Bag Optimizes Users Competencies in:

- Practicing handling and setting up of IV bags
- Understanding and adhering to established protocols and guidelines for IV therapy.
- Learning the correct procedures for hanging, priming, and adjusting the flow rate of IV bags.
- Developing safe handling techniques for IV bags, minimizing the risk of errors or contamination.
- Simulating patient interactions related to IV administration, fostering communication and bedside manner in a simulated environment
- Enhance critical thinking skills by making decisions related to IV therapy based on simulated scenarios.
- Simulating responses to potential emergencies or complications related to IV administration.

Qty: 1 bag

Product Content: Distilled water

NOT FOR HUMAN OR ANIMAL USE. FOR TRAINING PURPOSES ONLY.