



HELAGO-CZ, s.r.o.
Commercial Register maintained by the Regional Court in
Hradec Králové
Section C, File 17879
Kladská 1082
500 03 Hradec Králové 3
Company ID: 25 96 39 61, VAT: CZ 25963961
Phone: 495 220 229, 495 220 394
Fax: 495 220 154
GSM gate: 602 123 096
E-mail: info@helago-cz.cz
Web: http://www.helago-cz.cz

TFAW03 - Attachable Impaled Object simulator with bleeding set

Order code: **4134.TFAW03**



Cena bez DPH

301,00 Eur

Price with VAT

364,21 Eur

Parameters

Quantitative unit

ks

Introducing the Impaled Object Simulator, an innovative training device meticulously designed to replicate impaled object injuries on the body. This simulation tool aims to provide students and medical professionals with a realistic and immersive learning experience in

managing these critical and potentially life-threatening situations.

The Impaled Object Simulator features a two-part system: the first component simulates holed skin, while the second replicates a soft and flexible reinforcing rod. Together, these components create a lifelike representation of an impaled object scenario, enabling learners to practice their skills in a controlled and safe environment.

For optimal comfort and customization, the Impaled Object Simulator offers Our new system of attachment will help the users to adjust the simulator tight and comfortably on the body.

In summary, the Impaled Object Simulator stands as an exceptional training tool for honing skills in handling impaled object injuries. With its lifelike design, two attachment options, and focus on practicality, it provides a comprehensive and immersive training experience, empowering medical professionals and students to excel in critical situations and save lives.

Set Contents: Wearable wound, Squeeze bottle, Silicone tube, plastic connector and a plastic clamp

Size: L 28cm X H 0,6cm X W13,5cm

Weight: 170 grams